

# PATIENT-SPECIFIC SPECIALTY TOURNIQUETS

Successful venipuncture requires that a tourniquet be tied tightly enough to stop arterial blood flow, causing the collection vein to swell – so using a tourniquet with the right mechanical properties for a patient is absolutely necessary.

## Bariatric Tourniquets: Do Away with Tying Tourniquets Together

### THE CHALLENGE

When standard tourniquets are not long enough for larger patients, nurses will tie two tourniquets together to extend the length.

### THE SOLUTION

PlatinumCode bariatric tourniquets are twice as long as a standard tourniquet so that you can comfortably apply the band to the limb of a larger patient.

### THE RESULT

Bariatric tourniquets provide you with process efficiency while respectfully providing your patients with the comfort that they deserve.

### NO MORE TYING

Bariatric smooth tourniquets measure 1.5" x 32", twice the length and twice the width of standard tourniquets.



## Pediatric Tourniquets: Ready-to-Use Tourniquets for Toddlers and Children

### THE CHALLENGE

Standard tourniquets are often oversized for pediatric patients, so at many facilities today, nurses are cutting tourniquets to a shorter length.

### THE SOLUTION

PlatinumCode pediatric tourniquets are sized for little patients, allowing you to easily tie a tourniquet with the appropriate amount of pressure for venipuncture.

### THE RESULT

Pediatric tourniquets eliminate the need for cutting, giving you back time and reducing the risk of cross infection from a scissors.

### NO MORE CUTTING

Pediatric smooth tourniquets measure 0.75" x 12" and are specially designed to provide just enough grip and pressure for little ones.



## Infant Tourniquets: Use Products How They Were Made to Be Used

### THE CHALLENGE

NICU nurses, lacking an appropriately-sized solution for infant patients, use rubber bands or cut off glove tips for tourniquets.

### THE SOLUTION

PlatinumCode infant tourniquets are sized appropriately for infant limbs and come in accessible packaging and easily identifiable colors.

### THE RESULT

Infant tourniquets improve process efficiencies for NICU nurses, increase patient comfort, and reduce the risk of latex allergies and cross-contamination of non-tourniquet materials being used incorrectly.

### NO MORE RUBBER BANDS

Infant smooth tourniquets measure a wee 0.5" x 9" and come in unique "baby"-blue or purple colors.

