

Non-Latex Cohesive Bandages

Designed to stick to itself, rather than your patients' skin, Caresfield's cohesive bandages increase patient comfort and reduce skin sensitivity. With a wide variety of sizes and colors, including fun prints, you are sure to find your perfect fit!



PROVIDES PATIENT COMFORT AND SAFETY



- Bandages are self-adhering and don't stick to patients' skin
- Ideal for patients with sensitive skin
- Non-latex materials address risk of latex allergies
- Bright colors and pediatric prints promote a positive patient experience
- Single-use option provides just enough material for one patient, reducing risk of cross-infection

EASY TO USE



- Individually wrapped rolls
- Bandage stretch ensures proper coverage and pressure
- Wide variety of colors, sizes, and fun prints meet the needs of virtually every department in your facility



How To Tear Cohesive Bandages

Using a cohesive bandage is simple:

1

Pull Bandage



Pull desired amount of cohesive bandage out from roll. The cohesive bandage stretches and adheres to itself, therefore not much is needed.

2

Finger Placement



Place thumb and pointer finger close together at the very edge of the bandage as if to pinch the bandage. Quickly and with force, tear the cohesive bandage in opposite directions.

Standard Non-Latex Cohesive Bandages

Standard Colors								
Texture	Assorted	Tan	Pink	Orange	Green	Blue	Red	Sterile Tan*
1" x 5 yards	✓	✓	✓	✓	✓	✓	✓	✓
1.5" x 5 yards	✓	✓	✓	✓	✓	✓	✓	✓
2" x 5 yards	✓	✓	✓	✓	✓	✓	✓	✓
3" x 5 yards	✓	✓	✓	✓	✓			✓
4" x 5 yards	✓	✓	✓	✓	✓			
6" x 5 yards	✓	✓	✓	✓				
2" x 1 yards (Single-Use)	✓							

*Only the item numbers within this column are sterile.

Pediatric Non-Latex Cohesive Bandages

Fun Prints			
Texture	Camouflage	Pink Camouflage	Assorted Fun Prints
2" x 5 yards	✓	✓	✓
1.5" x 5 yards	✓	✓	
1" x 5 yards	✓	✓	

