

Leak-Resistant Specimen Transport Cups

Caresfield Specimen Transport Cups provide a safe, compliant, leak-resistant method of transporting specimen samples, whether that be across the hall or across the country.



Easy-to-read color changing temperature strip at cup base

Available in 90ml (3oz) or 120ml (4oz) cup volume

Sterile or non-sterile options (3 oz)



SMART FEATURES FOR LEAK RESISTANCE

- Click Top screw cap lets you know the cup is firmly sealed
- Tamper evident seals
- Leak-resistant lid rated to 95kPa for transport inside of pneumatic systems or via air (certified compliance with 49 CFR 178.605 and IATA DGR 6.3.5)



USER FRIENDLY

- Stackable design for storage space efficiency
- Translucent cup with graduated markings (ml and oz) for easy sample volume measurement
- Pre-applied I.D. label with name, D.O.B., date, and medical record number
- Ergonomic lid and base with ribbed closure facilitate easy gripping with gloves

Let's talk



952.890.7100



caresfield.com



sales@caresfield.com

Intuitive Design

When Opening



The lid seal will tear when opened to indicate sterility has been breached.

When Closing



Screw the lid on and listen for the "click" sound to ensure the cup is fully sealed.

When Testing



Look at the easy-to-read indicator strip to determine specimen temperature.

Click-Top Specimen Transport Cups

	ITEM NO.	DESCRIPTION	QTY.
	61405	90ml (3oz) Specimen Cup, Sterile, Tamper Label, Stackable	400/Case 100-cup packs, 4 packs per case
	61406	90ml (3oz) Specimen Cup, Sterile, Tamper Label, Temp Strip, Stackable	400/Case 100-cup packs, 4 packs per case
	61408	120ml (4oz) Specimen Cup, Sterile, Tamper Label, Stackable	400/Case 100-cup packs, 4 packs per case
	61409	120ml (4oz) Specimen Cup, Sterile, Tamper Label, Temp Strip, Stackable	400/Case 100-cup packs, 4 packs per case
	61239	90ml (3oz) Non-Sterile Urine Specimen Cup, White Cap	400/Case 100-cup packs, 4 packs per case

Specialty Specimen Cups



CLIA Screen In-Vitro Drug Test Cups are also available from Caresfield. For additional details, pricing, and availability, please contact us with your panel requirements by calling 952.890.7100 or visiting caresfield.com.